

## AVCPT FORMULA GUIDE

HEMATOLOGY		
Test	Formula	Example
Mean Corpuscular Volume (MCV)	$\frac{PCV \times 10}{RBC} = MCV \text{ (fL)}$	<p><b>Patient Information:</b>  <b>PCV = 45%</b>  <b>RBC = 7.50 x 10<sup>6</sup>/μL</b></p> <hr/> $\frac{45 \times 10}{7.50} = 60.0 \text{ fL}$ <p style="text-align: center;">MCV = 60.0 fL</p>
Mean Corpuscular Hemoglobin Concentration (MCHC)	$\frac{Hb \text{ (g/dL)} \times 100}{PCV} = MCHC \text{ (g/dL)}$	<p><b>PCV = 45%</b>  <b>Hb = 15 g/dL</b></p> <hr/> $\frac{15 \times 100}{45} = 33.3 \text{ g/dL}$ <p style="text-align: center;">MCHC = 33.3 g/dL</p>
Mean Corpuscular Hemoglobin (MCH)	$\frac{Hb \text{ (g/dL)} \times 10}{RBC} = MCH \text{ (pg)}$	<p><b>Hb = 15 g/dL</b>  <b>RBC = 7.50 x 10<sup>6</sup>/μL</b></p> <hr/> $\frac{15 \times 10}{7.5} = 20.0 \text{ pg}$ <p style="text-align: center;">MCH = 20.0 pg</p>
Red Blood Cell Estimate	$\frac{PCV}{6} = \sim \# RBC \times 10^6/\mu L$	<p><b>PCV = 45%</b></p> <hr/> $\frac{45}{6} = 7.5 \times 10^6/\mu L$ <p style="text-align: center;">RBC Estimate = 7.50 x 10<sup>6</sup>/μL</p>

Hemoglobin Estimate	$\frac{PCV}{3} = \sim \text{Hb g/L}$	<p><b>PCV = 45%</b></p> <hr/> $\frac{45}{3} = 15 \text{ g/dL}$ <p>Hb Estimate = 15 g/dL</p>
PCV Estimate	$\text{Hb} \times 3 = \sim \text{PCV (\%)}$	<p><b>Hb = 15 g/dL</b></p> <hr/> $15 \times 3 = 45\%$ <p>PCV Estimate = 45%</p>
Hemocytometer (Neubauer) Formula	$\frac{\text{Number of Cells Counted}}{\text{Number Large Squares Counted}} \times \text{depth of hemacytometer (10)} \times \text{dilution} = \text{cells}/\mu\text{L or mm}^3$	
Reticulocyte Percentage	$\frac{\text{Number Reticulocytes Counted}}{1000 \text{ RBCs Counted}} \times 100 = \% \text{ Reticulocytes}$	<p><b>Number of Retics / 1000 RBCs = 25</b></p> <hr/> $\frac{25}{1000} \times 100 = 2.5\%$ <p>Reticulocyte Percentage = 2.5%</p>
Absolute Reticulocyte Count	$\% \text{ Retics} \times \text{Total RBC count}/\mu\text{L} = \text{Retics}/\mu\text{L}$	<p><b>% Reticulocytes = 2.5%</b>  <b>Total RBC Count = 7,500,000/μL</b></p> <hr/> $\frac{2.5}{100} \times 7,500,000/\mu\text{L} = 187,000 \text{ Retic}/\mu\text{L}$
Corrected Reticulocyte Percentage (CRP)	$\% \text{ Retics} \times \frac{\text{Patient's PCV}}{\text{Avg. PCV for species}} = \% \text{ Reticulocytes}$ <p>Common "Normal" Avg. PCV: Canine = 45%</p>	<p><b>Patient Canine:</b>  <b>2.5% Retics</b>  <b>HCT: 25%</b></p> <hr/> $2.5 \times \frac{25}{45} = 1.39\%$ <p>CRP = 1.4%</p>

<p>Red Cell Production Index (RPI)</p> <p>Canine only</p>	$\frac{\text{CRP \%}}{\text{Maturation Factor}} = \text{RPI \%}$ <table border="1" data-bbox="720 277 1199 565"> <thead> <tr> <th>Patient PCV %</th> <th>Maturation Time (days)</th> </tr> </thead> <tbody> <tr> <td>45</td> <td>1</td> </tr> <tr> <td>35</td> <td>1.5</td> </tr> <tr> <td>25</td> <td>2.0</td> </tr> <tr> <td>15</td> <td>2.5</td> </tr> </tbody> </table>	Patient PCV %	Maturation Time (days)	45	1	35	1.5	25	2.0	15	2.5	<p><b>CRP = 1.4%</b> <b>HCT = 25%</b></p> <hr/> $\frac{1.4\%}{2.0} = 0.7\%$ <p>RPI = 0.7%</p>
Patient PCV %	Maturation Time (days)											
45	1											
35	1.5											
25	2.0											
15	2.5											
<p>Absolute WBC Counts</p>	<p>Total CWBC count/<math>\mu\text{L}</math> x (% of each WBC type/100) = absolute count of each cell type/<math>\mu\text{L}</math> or <math>\text{mm}^3</math></p>	<p><b>WBC count = 9000/<math>\mu\text{L}</math></b> <b>Neutrophils = 60%</b> <b>Lymphocytes = 35%</b> <b>Eosinophils = 5%</b></p> <hr/> $9,000 \times \frac{60}{100} = 5,400/\mu\text{L}$ <p>OR</p> $9,000 \times 0.60 = 5,400 \text{ neutrophils}/\mu\text{L}$ $9,000 \times \frac{35}{100} = 3,150/\mu\text{L}$ <p>OR</p> $9,000 \times 0.35 = 3,150 \text{ lymphocytes}/\mu\text{L}$ $9,000 \times \frac{5}{100} = 450/\mu\text{L}$ <p>OR</p> $9,000 \times 0.05 = 450 \text{ eosinophils}/\mu\text{L}$										

White Blood Cell Estimate (from peripheral blood film)  (High Dry Fields [10x ocular x 40 x obj. = 400x mag.]	$\frac{\text{number of WBC in 10 fields}}{10} \times 2,000 = \sim \text{number of WBC}/\mu\text{L}$	<b>50 WBCs in 10 HPFs (in the monolayer on 40x)</b>  $\frac{50}{10} \times 2,000 = 10,000/\mu\text{L}$
Corrected White Blood Cell Count (for the presence of Nucleated Red Blood Cells)	$\frac{\text{WBC count}/\mu\text{L} \times 100}{100 + \text{number of NRBC per 100 WBC}} = \text{corrected WBC count}/\mu\text{L}$	<b>WBC Count = 45,000/<math>\mu\text{L}</math> 20 NRBC per 100 WBC (from manual differential)</b>  $\frac{45,000 \times 100}{100 + 20} = 37,500/\mu\text{L}$
Platelet Estimate (from peripheral blood film)  (oil immersion fields [10x ocular x 100x obj. = 1000x mag.]	$\frac{\text{number of platelets in 10 fields}}{10} \times 20,000 = \sim \text{number of platelets}/\mu\text{L}$	<b>Platelets = 120/10 fields (in the monolayer on 100x)</b>  $\frac{120}{10} \times 20,000 = \frac{2,400,000}{10} = 240,000/\mu\text{L}$
Fibrinogen (Heat Precipitation Method)	(Non-incubated Plasma Protein g/dL minus incubated Plasma Protein g/dL) x 1000 = Fibrinogen mg/dL	<b>Non-incubated Plasma Protein = 6.0 g/dL Incubated Plasma Protein = 5.4 g/dL</b>  $(6.0 - 5.4) \times 1,000 = 600 \text{ mg/dL}$

## URINALYSIS

Test	Formula	Example
Dilution of USG when initial USG reading exceeds upper limits of the refractometer scale	Prepare a 1:2 dilution by mixing equal parts of urine and deionized / distilled water.  Read USG using a clinical refractometer. Multiply the digits to the right of the decimal by 2 to compensate for the 1:2 dilution.	<b>Initial reading exceeds the upper limits of the refractometer scale.</b>  USG on 1:2 diluted sample = 1.036 $0.036 \times 2 = 0.072$ $= \text{USG } 1.072$

## PARASITOLOGY

Test	Formula	Example
Fecal Egg Count Reduction Test (FECRT)	$\frac{\text{EPG (pre-treatment)} - \text{EPG (14 day post-treatment)}}{\text{EPG (pre-treatment)}} \times 100 = \%$	Pre-treatment on an equine = 1,050 epg Post-treatment (14 days) = 150 epg  $\frac{1,050 - 150}{1,050} = 86\%$
Modified McMasters Quantitative Fecal Exam (Eggs/g of feces [e/g])	Common Method of Ruminants with Large Volumes of Feces (example: cattle) using 4 g feces, 56 mL flotation Total number of eggs for each parasite counted x 50 = epg  OR  4 g feces, 26 mL flotation solution Total number of eggs for each parasite counted x 25 = epg	Patient cow: 4 g feces, 56 mL flotation Solution, total of 10 eggs counted in both chambers  10 x 50 = 500 epg
	Common Method for Ruminants with Small Fecal Volumes (example: sheep, goats) using 2 g feces, 28 mL flotation Total number of eggs for each parasite counted x 50 = epg	Patient sheep: 2 g feces, 28 mL flotation Solution, total of 5 eggs counted in both chambers  5 x 50 = 250 epg
	Common Method for Horses using 4 g feces, 26 mL flotation Total number of eggs for each parasite counted x 25 = epg	Patient horse: 4 g feces, 26 mL flotation Solution, total of 15 eggs counted in both chambers  15 x 25 = 375 epg

## REFERENCES

Cowell RL. *Veterinary Clinical Pathology Secrets*. Elsevier-Mosby, 2004: 11, 33, 35.

Latimar KS. *Duncan and Prasse's Veterinary Laboratory Medicine: Clinical Pathology*. Wiley-Blackwell, 5<sup>th</sup> edition, 2011: 13-15, 22-23, 59, 113.

Nielsen MK, et al. *AAEP Parasite Control Guidelines*. 2013; 4. ([www.aaep.org/custdocs/parasitecontrolguidelinesfinal.pdf](http://www.aaep.org/custdocs/parasitecontrolguidelinesfinal.pdf))

Rosenfeld AJ, Dial SM. *Clinical Pathology for the Veterinary Team*. Wiley, 2011; 28-29, 54-55.

Sink CA, et al. *Laboratory Urinalysis and Hematology for Small Animal Practitioner*. Teton New Media, 2004; 55-56, 62, 80-82, 98-99.

Sink CA, et al. *Practical Veterinary Urinalysis*. Wiley-Blackwell, 2012; 23.

Voigt GL, et al. *Hematology Techniques and Concepts for Veterinary Technicians*. Wiley-Blackwell, 2<sup>nd</sup> edition, 2011; 46-47, 59-60, 69, 117-119, 126-129.

Zajac AM, Conboy GA. *Veterinary Clinical Parasitology*. Wiley-Blackwell, 8<sup>th</sup> edition, 2012; 9-11.